

(800) 728-0974

CopQuest.com/nik

CopQuest™

ODV / NIK Crossover Chart

Presumptive Test	 Part No.	nik [®] Part No.	Test	CopQuest [®] Part No.
	Marquis Reagent — Opium alkaloids, Heroin and Amphetamine type compounds; general screening agent for other drugs	902	800-6071	A
Nitric Acid Reagent — Secondary screening test - Opiates (Morphine, Codeine or Heroin) and Amphetamine-type compound	903	800-6072	B	43-2101-000
Modified Dille-Koppanyi Reagent — Barbiturates; use after no change result in Test A and Test G, and orange result in Test I	905	800-6073	C	43-2102-000
Modified Ehrlich's Reagent — LSD (D-Lysergic Acid Diethylamide); use after positive results in Test A and B	907	800-6074	D	43-2103-000
Duquenois-Levine Reagent — Standalone test for Marijuana, Hashish and Hash Oil	908	800-6075	E	43-2104-000
Acid Neutralizer — Use upon completion of each test	910	800-6076	F	43-2105-000
Modified Scott Reagent — Cocaine, Crack or Free Base	904B	800-6077	G	43-2106-000
Methadone	906	800-6078	H	43-2107-000
Lieberman's Reagent — Screening of PMA, Ketamine, Barbiturates and Methadone	N/A	800-6089	I	43-2108-000
PCP (Phencyclidine). Test J used after no change results from Tests A, G, I and W	914	800-6079	J	43-2109-000
Heroin, Black Tar, Codeine and Morphine	922	800-6080	K	43-2110-000
Modified Mecke's Reagent — Heroin, including White, Brown and Black Tar, and MDMA (Ecstasy)	924	800-6081	L	43-2111-000
Methaqualone (Quaaludes, Sopor Somnifac, Opitamil and Parest)	914	800-6082	M	43-2112-000
Pentazocine (Talwin Nx or Talacen); Oxycodone	926	800-6083	N	43-2113-000
GHB (Gamma-Hydroxybutyrate)	928	800-6090	O	43-2114-000
Propoxyphene (Darvon)	N/A	800-6084	P	43-2115-000
Ephedrine and Pseudoephedrine	927	800-6085	Q	43-2116-000
Valium (Diazepam), Rohypnol (Flunitrazepam) and Methcathinone	925	800-6086	R	43-2117-000
Marijuana, green leaves, seeds	909	T/C	S	43-2121-000
Morris Reagent — Ketamine	912	800-6091	T	43-2120-000
Methamphetamine and MDMA (Ecstasy)	923	800-6087	U	43-2118-000
Mandelin Reagent — Amphetamines and Methadone; screening for PMA and Ketamine in conjunction with Test I	906	800-6088	W	43-2119-000
Mayer's Reagent	901	N/A	N/A	(None)

Information here is from the manufacture. It is believed to be accurate at the time of publishing but is not guaranteed.